


---

**PATIENT INFORMATION  
&  
INSTRUCTIONS**

---

**The "Millie"™ Universal Urinal**



PATENTED - U.S. PAT #4665571

**For Men and Women**

For more information contact:

**VLSCOT**™  
MEDICAL LLC

P.O. Box 351, East Hanover, NJ 07936

Tel. (973) 887-9273

Fax (973) 887-3961

**Outside NJ (800) 221-0658**

[www.viscot.com](http://www.viscot.com)

For more information contact:

**VLSCOT**™  
MEDICAL LLC

P.O. Box 351, East Hanover, NJ 07936

Tel. (973) 887-9273


Fax (973) 887-3961

**Outside NJ (800) 221-0658**

[www.viscot.com](http://www.viscot.com)


## INSTRUCTIONS

1. Grasp the Urinal by the neck with the handle side up and the opening toward the body.



(**Important:** Handle must be up for the Urinal to fit against the body properly and to avoid spillage)

The lip of the Urinal has been specially designed to prevent spillage.



2. Bring the opening of the Urinal up firmly against the body.



3. When finished using the Urinal, simply hook the handle onto the bed rail or chair for measurement, disposal and easy access.

**Note:** This Urinal is intended to be used by only one patient.

**Note:** This Urinal may be disinfected with a 50% solution of bleach and water then easily rinsed with water.

While the "Millie" Universal Urinal can be used by everyone, it is uniquely designed to fit the needs of the female patient. When used properly, it can aid the female patient to urinate in a sitting or lying down position while in bed without spillage. The patient is easily able to use the "Millie" Universal Urinal conveniently and in complete privacy, without assistance.

The neck of the "Millie" Universal Urinal is designed to be easily grasped with one hand. The opening is tapered to be anatomically correct.

The easy to use and comfortable "Millie" Universal Urinal has a graduated scale for measurement and visual inspection. Also, a convenient handrail for hanging the unit on the bed rail or a chair. The urinal can be easily transported for disposal of contents without spillage.