

# CHG Compatible *VS* Prep Resistant



Withstands up to a 30 second prep with CHG prep solutions.



Increased Gentian Violet Dye used in formula.



Contains alcohol and can be removed with alcohol based preps.



Only available in purple.



Clearly visible after 60 seconds of prep with CHG preps.



Unique patented formula engineered to stand up to *all* common prep solutions.



Alcohol free to stand up against alcohol based prep solutions.



Available in purple, green, and red to maximize visibility on all skin tones.

Not all "prep resistant" inks are the same.  
Let us help you see the difference.